



The Advantages at a Glance

- quick fabrication of bases and individual trays
- exact fitting after polymerization
- insulation of plaster models is not required

Range of Use

- total and partial prosthetic technique
- individual trays
- bases for bite patterns and odontograms

Processing Instructions

- Mould the light-curing base plates onto the model and cut off the excess material.
- Form the material labial/buccal to the jaw and smooth it by hand. In case of partial models cut out the base dependent on the teeth.
- The better the modelling is done, the less time for proceeding is required.
- Polymerize the plates for around 5 minutes under UVA-light and then finish.
- To get back the original shade, apply vaseline or monomer on whitish shining grinding areas.
- Store cool (<20°C/68°F) and light protected.

Presentations and Prices

Item	Number	Item N°	€/50 pcs.	Number	Item N°	€/100 pcs.
pink upper	_____	90101	41,-	_____	90102	79,-
pink lower	_____	90103	41,-	_____	90104	79,-
blue upper	_____	90105	41,-	_____	90106	79,-
blue lower	_____	90107	41,-	_____	90108	79,-
transparent- upper	_____	90109	43,-	_____	90110	83,-
mintflavoured lower	_____	90111	43,-	_____	90112	83,-
natural upper	_____	90113	41,-	_____	90114	79,-
natural lower	_____	90115	41,-	_____	90116	79,-
green- upper	_____	90117	43,-	_____	90118	83,-
mintflavoured lower	_____	90119	43,-	_____	90120	83,-

Fax-Order

Fill in required numbers above and send to (+ +49) (0)234-33 0545.
For multiple orders please copy this form.

Order placed by: _____

address: _____

phone: _____

fax: _____

date: _____

stamp and signature: _____

Prices net of VAT, packing and transport.

Technical alterations reserved.

Delivery is subject to our General Terms and Conditions.



Clemensstraße 14
D 44 789 Bochum
fon: (+ +49) (0)234-30 17 84
fax: (+ +49) (0)2 34-33 05 45
web: www.vigano-dental.de
email: info@vigano-dental.de